



University Hospital
Southampton
NHS Foundation Trust

Review of expected results for mA variation curves in routine CT AEC testing

CT Users Group Annual Meeting

16th October 2025

Oliver Lawes
Trainee Clinical Scientist

Contents

- Background & Method
- GE & Siemens Results
- Lungman Results
- Canon Results
- A question to the audience...
- Quantitative assessment
- Limitations & Conclusions


Background

- CT AEC testing is performed using the nested CTDI phantom method at UHS [1]
- AEC reproducibility is confirmed by graphing the mA variation along the z-axis and visually assessing the deviation from baseline
 - Performed only to confirm the AEC is functioning, so has no suspension criteria (performance assessed by dose audits)
- It was noticed that results from GE scanners tended to be infrequently reproducible
 - Was this operator error? Was the method unsuitable? Was the scanner incapable of reproducible results?

[1] Iball, G.R., Moore, A.C. and Crawford, E.J. (2016), A routine quality assurance test for CT automatic exposure control systems. Journal of Applied Clinical Medical Physics, 17: 291-306. <https://doi.org/10.1120/jacmp.v17i4.6165>

Background

- CT AEC test at UHS [1]
- AEC reproduction along the z-axis
 - Performance criteria
- It was noticed that the test was not reproducible
 - Was this due to the scanner's automatic exposure control system being infrequently calibrated?



phantom method

variation along
baseline

no suspension

be infrequently



s the scanner

[1] Iball, G.R., Moore, A.C. and Crawford, E.J. (2016), A routine quality assurance test for CT automatic exposure control systems. Journal of Applied Clinical Medical Physics, 17: 291-306. <https://doi.org/10.1120/jacmp.v17i4.6165>

- Repeatability \approx Reproducibility
 - Perform the standard AEC test 3 times back to back and compare results
 - Repeat on a Siemens scanner to compare with GE
- Would an anatomical phantom produce more consistent results?
 - Perform the standard AEC test but using a Kyoto Kagaku Lungman phantom instead of the CTDI phantom
 - Chest plates were used





Results – GE



Protocols - GE

	AEC Baseline Protocol Set 1 Fill out Baselines tab with all planned protocol information and will auto-populate below.
Protocol name	7.1 Pelvis Prosoma
Pitch	1.375
AEC setting	NI = 26.45
kV	120
mA range	100 - 440
Rotation time (s)	1
Collimation	1.25 x 16i
Slice width (mm)	1.25
FOV	large (50)
Filter (H/B)	body
Scan direction	Head first
Table height	261
Reconstruction type (i.e recon kernel, iterative, % etc)	std, full, 50% ASIR



Results – Siemens



Protocols - Siemens



	AEC Baseline Protocol Set 1 Fill out Baselines tab with all planned protocol information and will auto-populate below.
Protocol name	AbdomenRoutine
Pitch	0.6
AEC setting	CAREdose4D on, CAREkV off, Ref mAs=147, dose opt for 7
kV	120
mA range	17-555
Rotation time (s)	0.5
Collimation	128x0.6
Slice width (mm)	5
FOV	400mm
Filter (H/B)	Standard
Scan direction	Craniocaudal
Table height	159
Reconstruction type (i.e recon kernel, iterative, % etc)	Br38

Results – Lungman





- No improvement in repeatability over the CTDI
- Offered no benefit to visual assessment either

Results – Canon



Results – Canon




Protocols - Canon

	AEC Baseline Protocol Set 1 Fill out Baselines tab with all planned protocol information and will auto-populate below.
Protocol name	User Abdo Non-contrast
Pitch	0.813
AEC setting	SureEXP 3D, sd=10
kV	120
mA range	80-500
Rotation time (s)	1
Collimation	0.5 x 80
Slice width (mm)	0.5
FOV	320 (M)
Filter (H/B)	M
Scan direction	Craniocaudal
Table height	122
Reconstruction type (i.e recon kernel, iterative, % etc)	XY on, OE off

Time for a vote

- Would you say that the Siemens results were the most reproducible?





Results – Quantitative

- Maximum variation between 2 individual points from:
 - Anywhere on the graphs
 - Only from points across a plateau

Max Variation	Reproducibility whole curve	Reproducibility plateau	Repeatability whole curve	Repeatability plateau
GE (2021, 2023)	24%	9%	62%	11%
Siemens	42%	8%	29%	5%
Canon (excl. 2021)	23%	12%	5%	5%

- The aim was to derive a tolerance to replace visual check
 - An average variation across whole curve would lead to max variations being out of tolerance

Results – Quantitative



Results – Quantitative

- Percentage tolerance based on the overall reproducibility deviation is too large
 - Based on the deviation across plateaus is too small
- What is a suitable value for a fixed tolerance to encompass entire graph?
- Excel isn't suitable for graphing a tolerance
 - Starts graph at incorrect z position
 - Doesn't work for non integer slice divisions

Results – Quantitative

- What about an R^2 tolerance between the curve and baseline?
 - Less of a visual aid, more of a confirmation

GE

RSQ	2020	2021	2022	2023	2024
2020	1.00				
2021	0.78	1.00			
2022	0.00	0.04	1.00		
2023	0.86	0.98	0.02	1.00	
2024	0.21	0.42	0.64	0.37	1.00

RSQ	1	2	3
1	1.00		
2	0.88	1.00	
3	0.86	0.99	1.00

Siemens

RSQ	2020	2021	2022	2023
2020	1.00			
2021	0.97	1.00		
2022	0.90	0.97	1.00	
2023	0.99	0.99	0.93	1.00


RSQ	1	2	3
1	1.00		
2	1.00	1.00	
3	0.99	0.99	1.00

Canon


RSQ	2021	2022	2023	2024
2021	1.00			
2022	0.89	1.00		
2023	0.84	0.85	1.00	
2024	0.81	0.84	0.72	1.00

RSQ	1	2	3
1	1.00		
2	1.00	1.00	
3	1.00	1.00	1.00

Results – Quantitative




RSQ	2020	2021	2022	2023	2024
2020	1.00				
2021	0.78	1.00			
2022	0.00	0.04	1.00		
2023	0.86	0.98	0.02	1.00	
2024	0.21	0.42	0.64	0.37	1.00




RSQ	1	2	3
1	1.00		
2	0.88	1.00	
3	0.86	0.99	1.00

Results – Quantitative




RSQ	2020	2021	2022	2023
2020	1.00			
2021	0.97	1.00		
2022	0.90	0.97	1.00	
2023	0.99	0.99	0.93	1.00




RSQ	1	2	3
1	1.00		
2	1.00	1.00	
3	0.99	0.99	1.00

Results – Quantitative



RSQ	2021	2022	2023	2024
2021	1.00			
2022	0.89	1.00		
2023	0.84	0.85	1.00	
2024	0.81	0.84	0.72	1.00

*2020 omitted as less data points



RSQ	1	2	3
1	1.00		
2	1.00	1.00	
3	1.00	1.00	1.00

Results – Quantitative

- R^2 comparison does not aid the visual check
- Nor can it produce a reliable tolerance from repeatability results to provide a confirmation of visual judgement
 - Reproducibility results could potentially be used for a system with reliable reproducibility, but unhelpful at acceptance / without historical data

Conclusions

- Visual assessment is sufficient and simpler than quantitative
- Lungman offers no benefit over the CTDI phantom
- Key criteria can be identified for GE and Siemens:
 - GE: Overall shape should be consistent, mA variations are acceptable across some but not all sections
 - Siemens: Overall shape and mA should be consistent
- Canon requires more work to determine similar guidance criteria
- *We aren't certain on our method however, is anyone doing the same or anything differently?*

Limitations

- Only one scanner of each manufacturer used
- Only one protocol (our baseline protocol) used
- Only mA considered, not noise or CT number
- Did not vary protocol setup to try and identify reasons for historical deviations
 - Suspect larger variations in GE results is from user mistakes
- Have not confirmed whether new / future results match the identified criteria



University Hospital Southampton
NHS Foundation Trust

Thank you for listening

Questions?

University Hospital Southampton
NHS Foundation Trust
Tremona Road, Southampton
Hampshire, SO16 6YD

www.uhs.nhs.uk

Oliver.Lawes@uhs.nhs.uk